Male External Catheter Size Guide

How to determine your Male External Catheter Size

Print this PDF – When printing, please ensure any page scaling is turned off by selecting Print Actual Size (or similar depending on your print program). Once printed, check image is correct size by aligning \$1 coin to image shown.

Option 1:

Cut out the guide located on the right side of this page along the dotted line. Place the guide over the flaccid penis behind the head and choose the size that fits best.

Option 2:

Cut out the measuring tape, and place it around the flaccid penis behind the head. Measure the circumference of the penis and then divide the result by 3.14. Eg: 12cm / 3.14 = 3.8cm (38mm).


If you are still not sure which size will suit, we suggest you order a Paralogic sample pack which includes 3 catheters in different sizes to determine your best fit - go to:

http://www.paralogic.com.au/male-external-catheters

If you still have doubts we recommend you order the smaller size.

The catheters are made of silicone and have some stretch and it is better to avoid leaks.









1300 727 222 www.paralogic.com.au orders@paralogic.com.au